



# Healthcare Accessories

## Patient Chair/Recliner Options

### Folding Side Table

Although the optional Folding Side Table must be ordered factory installed, however, some assembly is required. Remove shelf with brackets from the inner carton. Locate mounting holes on the outside of the chair/recliner arm. Align the bracket with the mounting holes in the arm and insert bolts from hardware packet. Be sure to tighten completely.



### Pivoting Tablet

Although the optional Pivoting Tablet must be ordered factory installed, some minor assembly and adjustment is required. Remove the tablet arm and Allen wrench (for tensioning adjustment) from the inner carton. With the (enclosed) Allen wrench, loosen the tensioning screw on the side of the mount by turning counter-clockwise, but do not remove. Place the black stem into the mount on the side of the unit's arm. Once the tablet arm is fully seated into the mount, tighten the set screw (clockwise) until the desired tension is experienced. Note: If tightened too much, the arm will not rotate.

### Push Bar

Although the optional Push Bar must be ordered as factory installed, some assembly is required. Locate the three mounting holes on each side of the upper back of unit. Carefully place the push bar in position aligning the holes in the push bar bracket with the holes in the upper back. Start one bolt on each side, but do not tighten yet. Now start the other two bolts on each side of back. Now tighten all six completely using the Allen wrenches provided in pack.



### IV Pole

Although the optional IV Pole must be ordered as factory installed, some assembly is required. Remove IV pole and hanger bar from shipping carton. Tighten hanger bar onto top of IV pole. Place assembled pole into mount on bottom rear of the arm. Tighten locking knob to secure IV pole.