



Hannah Treatment Chair



Goes Great with Work.®



Comfort, form and function are the essential building blocks for any HPFI® product. Nowhere is this more critical than in a healthcare setting where patient comfort is essential to maximize relaxation and recovery while simultaneously providing caregivers with equipment that looks great and performs flawlessly.

Recovery has never felt so good.

Standard Features

- Lay-flat mechanism allows for administering a wide variety of treatments without use of a bed or gurney
- Independently operated back and foot rest
- Pneumatically controlled back with infinite adjustment
- Sinuous wire spring suspension and high density foam for durability and patient comfort
- Full length foot rest provides support to patients entire leg
- Pull-out foot tray for easy patient transfer*
- Swivel casters (front) and locking swivel casters (rear) offer zero-turning radius for maximum maneuverability
- Back release located on outside of arm, right facing standard
- Standard width models rated to 500 lbs.
- Pillow back models include repositionable pillow.
- Warranty: Limited lifetime on seating components. Three years on wearable parts (recliner mechanism and casters)



Trendelenburg Tilt Feature

Offers a negative 3 degree pitch from the full flatbed position*

*Trendelenburg feature and pull out foot tray included on Models 830, 832, 834, 836

Additional Options

- IV pole and holder
- Back release can be located on inside of either arm or outside left facing
- Two color upholstery combinations available
 - Seat/Back: First fabric
 - Arms: Second fabric
 - (Cover photo fabric combination is a custom quote)
- Clear vinyl headrest covers for non-pillow back models

Dimensions

- 33"W (Models 830, 831, 834, 835)
- 37"W (Models 832, 833, 836, 837)
- 39"D (73" Reclined flat)
- 44"H



Brushed Aluminum or Black Powder Coated Push Bar

Optional black or brushed aluminum push bar with black textured grip.



Chair Controls

Optional brushed aluminum foot rest lever with black textured grip to replace standard black lever. Hands free foot pedal to activate Trendelenburg feature.



Tablet Arms

Select from two optional tablet arms.



Central Braking

Four heavy-duty central braking swivel casters controlled with three position pedal located on left side facing.



High-Carbon Steel Frame

Delrin bushings at critical wear points. Wear-tested at 50,000 human occupant cycles.



Tente Casters

3" non marking dual wheel swivel casters. Pressed steel housing with double ball bearing swivel head. Locking rear casters.



Pillow Back

Tufted look pillow back on 834, 835, 836 and 837.



Heat & Massage

Operated through push-button control located on inside of arm (not available on pillow back models).



Arm Caps

Black polyurethane or wood arm caps (available in 11 standard finishes.)



Folding Side Table

Optional side table locks solidly in place yet easily folds flat. Shown on Haley Recliner 840.



Bariatric Widths

28.5" inside clearance rated at 750 lb. capacity.



Accessory Hook

Provides convenient placement for drainage bags or securing restraints.

Shown on front:

830

Seat/Arms: Foundation 10 350-000 Chocolate

Back: Round N' Round WC762-007 Sandpiper CRYPTON®

